


stryker®

Craniomaxillofacial



Surgical Drains

**High-quality disposable operating room
products for wound drainage**

SPG Bulb Drainage System



For procedures requiring a small to mid-range reservoir, the SPG Bulb Drainage Kit provides convenient options for plastic, aesthetic and reconstructive surgical procedures.

The SPG Bulb Drainage System:

- Designed to provide gentle suction and prohibits reverse flow by relying on a 100cc silicone bulb reservoir with an anti-reflux valve for constant, steady evacuation of post-operative fluids
- May simplify ordering by providing conveniently packaged, high quality disposable products for wound drainage
- Packaged sterile
- Contain a radiopaque silicone catheter, trocar, and reservoir

TLS Surgical Drainage System



For surgeries of the craniofacial region, extremities, and other small fluid volume procedures.


The TLS Surgical Drainage System:

- Can be incorporated into the wound dressing
- Designed to provide convenient and effective closed suction drainage
- Utilizes standard, 15cc evacuated blood collection tubes as reservoirs
- Packaged sterile
- Contains a radiopaque silicone catheter, trocar, protective housing, connecting adapter, and two 15cc evacuated tubes

Kit contains: 1 Catheter, 1 Trocar, 2 Evacuated Tubes and a Protective Housing


CAT#	DESCRIPTION	
6640	7 FR	6/PK
6650	10 FR	6/PK
6660	4mm Flat with Side Holes	6/PK

TLS Surgical Drain Evacuation Tube Replacement Instructions




STEP 1

When evacuation tube is three-quarters full, slide clamp to closed position over drainage tubing.



STEP 2

Grip clear housing and white hub; remove evacuation tube and hub assembly from housing.




STEP 3

Grip cap of evacuated tube thru slots; then grip white hub and gently pull to remove. Be careful of needle.


Kit Contains: 1 Catheter, 1 Trocar, and 1 Reservoir

CAT#	DESCRIPTION
6646	15 FR, Round, 100cc Bulb 12/PK
6647	4mm Flat w/Side Holes, 100cc Bulb 12/PK
6648	10FR, Round, 100cc Bulb 12/PK
6649	4mm Flat w/Center Holes, 100cc Bulb 12/PK


All illustrations are shown at 400%.




4mm Flat with
center holes



4mm Flat with
side holes



10FR/3.33mm




15FR/5.0mm

QUINTUBE Monitor Pack




The QUINTUBE Monitor Pack, consisting of five evacuated tubes for purposes of monitoring drainage volume and color over time, is available as replacement reservoirs for the TLS Surgical Drainage System.

CAT#	DESCRIPTION
6630	QUINTUBE Monitor Pack, 5 Evacuated Tubes/Pack, 24 Packs/Case




STEP 4

Place new evacuation tube on white hub; push to puncture needle through center of tube stopper.



STEP 5

Replace hub and tube assembly into housing. Secure assembly near drainage location.



STEP 6

Slide Clamp to open position to re-institute drainage suction.



Joint Replacements

Trauma, Extremities & Deformities

Craniomaxillofacial

Spine

Biologics

Surgical Products

Neuro & ENT

Interventional Spine

Navigation

Endoscopy

Communications

Imaging

Patient Care & Handling Equipment

EMS Equipment

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: **QUINTUBE, SPG, Stryker, TLS**. All other trademarks are trademarks of their respective owners or holders.

Literature Number: **9410-400-216 Rev. A**
UnDe/P.S.

Copyright © 2011 Stryker
Printed in USA

Stryker Craniomaxillofacial
Kalamazoo, MI 49002 USA
t: 269 324 5346
toll free: 800 962 6558
f: 877 648 7114

www.stryker.com